


ACCEPTABLE FILL LEVELS FOR COMMON HAEMATOLOGY BOTTLES


STATED SAMPLE VOLUME

EDTA


1.3ML
RED CAP
(FBC)



3.0ML
MAUVE /
PINK CAP




1.3ML
GREEN CAP
(COAG)




CITRATE

2.7 / 4.5ML
BLUE CAP
(COAG)



5.0ML
BLACK
(ESR)



Fill level for coagulation samples is critical
Using Vacutainer system will fill bottles to required levels